Clay mineral analysis by XRD

From sample preparation to evidence review: explore the methodology guiding diffractogram comparison in Argiloteca.

This guide combines reading summaries with the full revised methodology. The described reference set uses standard samples from LDRX/CPGq/IGEO/UFRGS.

The PDF includes summaries and every full-text section, including collapsed sections. The original text remains in Portuguese.

Sample preparation · AI-generated conceptual illustration. It does not depict the laboratory or an analysed sample.

Research prototype · technical review required. The workbench organises curves and evidence to support interpretation. Similarity and scores do not represent mineral content, calibrated probability or independent validation.

On this page

XRD principles 2

N/G/H preparation and acquisition 4

Data integrity and organisation 7

Processing and peaks 8

Mineralogical interpretation 10

Curve similarity 16

Crystallographic references 17

Quantification and other techniques 19

Knowledge and artificial intelligence 20

Training, testing and validation 20

Using the workbench and preserving results 22

References and original document 23

Read each topic summary or expand “Full text” to consult the original wording and tables. Parameters describe the documented procedure; check each execution’s version and provenance.

Full text · original introduction in Portuguese

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XRD principles


Lamellar habit · AI-generated illustration without a scale. It is not a micrograph and does not identify a mineral species.

XRD records diffracted intensity against 2θ. Peak positions, intensities, shapes and widths must be interpreted with the preparation, radiation and processing used.

Bragg’s law relates angular position to interplanar spacing. In the v4 execution described in the document, a curve without an explicitly declared wavelength stays in 2θ, without conversion to d. The Cu Kα reference of approximately 1.5406 Å does not replace checking the actual radiation.

From instrument to diffractogram
An illustrative example of a diffractometer; the photograph does not depict the Argiloteca laboratory.
Sufangxi · Wikimedia Commons · CC BY-SA 4.0

The geometry of diffraction
The diagram relates the beam angle to the spacing between crystal planes.
Matthias M. / Dipl.-Phys. · Wikimedia Commons · Public domain

nλ = 2d sin θ

Bragg’s law · λ: wavelength; d: interplanar spacing; θ: half of 2θ; n: diffraction order. The described procedure uses n = 1.

Full text in Portuguese · XRD principles

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N/G/H preparation and acquisition

N is the natural preparation; G is glycolated; H is heated. H corresponds to the C (calcined) condition used in the workbench. N/G/H and N/G/C represent the same three conditions; original file labels are retained.

The protocol describes drying below 60 °C, disaggregation and quartering; the fraction below 2 µm is separated by sedimentation for oriented mounts. Unoriented preparation using a 250-mesh sieve is a separate procedure.

The suspension is sonicated for five minutes; after 2 h 30 min of settling, the upper 3 cm are collected. The stated centrifugation is 15,000 rpm for five minutes; recording the rotor or relative centrifugal force supports comparison between instruments. The paste is redispersed in 5–10 mL and dried on a glass slide for 24–48 hours.

Ethylene-glycol treatment takes place at 60 °C for two to three hours; H/C heating uses 550 °C for two hours. An oriented scan ending at 32° does not cover the 060 region with Cu Kα; this evidence requires the wider powder acquisition.

N · Natural

Reference condition before treatment.

G · Glycolated

Ethylene glycol · 60 °C · 2–3 h.

H/C · Heated

550 °C · 2 h.

Usual conditions described for the SIEMENS-BRUKER-AXS-D5000

Preparation	2θ	Step	Time per step
Unoriented powder	2°–72°	0,02°	1 s
Natural · N	2°–32°	0,02°	2 s

Preparation	2θ	Step	Time per step
Glycolated · G	2°–32°	0,02°	3 s
Heated · H/C	2°–32°	0,02°	2 s

Cu Kα radiation · 40 kV · 30 mA. Raw counts with different counting times require acquisition-time consideration; do not repeat conversion when the file already provides counts per second.

Full text in Portuguese · N/G/H preparation and acquisition

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Data integrity and organisation

Ingestion checks the signature and layout, header, point count, numeric values and increasing 2θ axis. The RAW extension alone does not guarantee support. Truncated files and inconsistent series require review.

The original is preserved and identified by SHA-256, format, size, angular range and reader version. Anonymised public copies use ARC codes; originals and identity mappings remain restricted. Recovering or classifying a curve does not authorise publication.

Full text in Portuguese · Data integrity and organisation

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Processing and peaks

The document describes ASLS background correction, Savitzky–Golay smoothing, robust noise estimation, peak detection and characterisation, and normalisation. Any alternative method must be recorded in the execution together with its parameters.

The ASLS smoothness parameter is not the radiation wavelength. Refining a maximum by parabolic interpolation does not correct instrument calibration. Local area uses trapezoidal integration; half-prominence width and background-referenced FWHM can differ.

Normalisation by the maximum, area or a region of interest facilitates shape comparison. None of these transformations provides mineral content. The parameters below are those documented; each result must retain its execution version and configuration.

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Mineralogical interpretation

The Å windows define search regions in the argiloteca.drx.ngc.v3 catalogue. They are project operating parameters, not universal limits from the literature. The semantic catalogue contains 18 rules, 34 windows, eight peak sets and six ambiguity groups; this does not establish uniform execution of every criterion or validation of 18 minerals.

Interpretation combines the N/G/H trajectory, reflection series and competing hypotheses. Smectites require investigation of glycolation expansion and thermal contraction; the 10 Å region may contain mica/illite and heated smectite. Kaolinite and chlorite may overlap near 7 Å. A single peak does not confirm a mineral.

The documentation distinguishes interstratification from a physical mixture, retains vermiculite as a provisional hypothesis and requires additional evidence to refine species. The 060 reflection requires adequate angular coverage and control of quartz interference near 1.54 Å. An unmeasured region is not equivalent to a reflection sought but not detected.

Catalogue operating windows · spacing in Å

Group or hypothesis	N	G	H/C
Smectites	13,46–16,86	16,06–18,31	9,65–10,37
Chlorites	13,74–14,74	13,83–14,72	13,58–14,87
Illite or mica	9,84–10,36	9,82–10,30	9,73–10,38

Group or hypothesis	N	G	H/C
Kaolinite group	6,97–7,42	6,96–7,42	Reduction or loss within 6.96–7.42
Vermiculites	13,70–14,80	13,70–14,80	10,00–12,00

For reference sample 7434, the recorded phases are smectite, mica/illite and kaolinite. The spreadsheet does not contain curves demonstrating treatment responses; those behaviours must be checked in linked files.

Full text in Portuguese · Mineralogical interpretation

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Curve similarity

The described similarity procedure restricts curves to their shared angular range and interpolates up to two thousand points without extrapolating unmeasured regions. It combines Pearson correlation (45%), an RMSE-derived measure (30%) and peak matching (25%).

Pattern similarity and the N/G/C evidence score are different results. Neither represents calibrated probability, purity or content. The compatible RAW1.01 administrative workflow uses the header axis without automatic shifting; other workflows must be interpreted through their own provenance.

Full text in Portuguese · Curve similarity

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Crystallographic references

A CIF describes a structural model. Patterns calculated with pymatgen are idealised references, distinct from measured diffractograms. The crystallographic pilot operates in shadow mode and adds candidates without replacing the primary classification.

The document records a local index of 12 structures for 10 minerals, calculated for λ = 1.5406 Å and 2θ from 2° to 70°. COD, the COD Inorganics subset and AMCSD are distinct sources; their full collections were not incorporated into the index. Public V1 uses peak positions; experimental V2 also considers intensities.

Full text in Portuguese · Crystallographic references

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Quantification and other techniques

Identification investigates compatible phases; quantification estimates proportions using preparation, models or standards, calibration and uncertainty analysis. Rietveld refinement, full-pattern fitting and internal standards are analytical options, not automatically validated quantitative capabilities of Argiloteca.

XRF, infrared spectroscopy, thermal analysis and electron microscopy may help discriminate hypotheses. Every experiment must remain associated with the same sample, its conditions and limitations.

Full text in Portuguese · Quantification and other techniques

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Knowledge and artificial intelligence

The vocabulary organises names and relations between samples, minerals, reflections, treatments, rules and interpretations. Terms such as “near”, “persistent” and “strongly reduced” require operating definitions before becoming computational criteria.

AI can support documentation, knowledge organisation and development. Scientific modules execute processing, while interpretive decisions remain human. A chat interface does not establish a working repository RAG connection or validated fine-tuning; evidence retrieval requires traceable sources.

Full text in Portuguese · Knowledge and artificial intelligence

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Training, testing and validation

The standard samples from LDRX/CPGq/IGEO/UFRGS form the reference set described for tests and training. Identifiers, curves, preparations, acquisition and mineral associations must be linked. A sample may contain multiple phases; Phase 1, Phase 2 and Phase 3 are neither treatments nor proportions.

The described spreadsheet has 54 entries and 53 distinct identifiers, including a repeated 7437. Identifier 106186 repeats smectite in its first two phases; X, ? and ? C require a legend. These values must not automatically become absence, zero or new classes.

Training, tuning and testing must be split by sample or source material, keeping preparations and derivatives together. Exploratory trials and software tests do not estimate production mineralogical accuracy. Independent validation requires experts, representative materials and analysis of errors and inconclusive outcomes.

Full text in Portuguese · Training, testing and validation

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Using the workbench and preserving results

Select analyses or upload temporary files, check the sample, preparation, radiation and angular range, compare curves, and examine evidence and interference. Expert review records the responsible person, date, analysis version and rationale; the conclusion may remain inconclusive.

Temporary files serve the session. Preparing a draft, exporting a report and publishing in InvenioRDM are separate operations. The N/G/C administrative package contains interactive HTML, JSON, CSV with original points and intensities, provenance and licence information; the HTML can be printed to PDF.

1. Select analyses or temporary files.
2. Check the sample, preparation, radiation and angular range.
3. Compare curves, peaks and competing hypotheses.
4. Record the review and preserve provenance.

Full text in Portuguese · Using the workbench and preserving results

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References and original document

The full revised document is available in Portuguese in the expandable sections and complete PDF document. References support the criteria; numerical windows and local weights belong to the Argiloteca procedure.

Bibliographic foundation

- Brindley, G. W.; Brown, G. (eds.), 1980. *Crystal Structures of Clay Minerals and their X-Ray Identification*. Mineralogical Society Monograph, 5.
- Bailey, S. W., 1980. *Structures of Layer Silicates*. Chapter in *Crystal Structures of Clay Minerals and their X-Ray Identification*.
- Bailey, S. W. (ed.), 1988. *Hydrous Phyllosilicates (Exclusive of Micas)*. Reviews in Mineralogy, 19.
- Moore, D. M.; Reynolds, R. C., Jr., 1989/1997. *X-Ray Diffraction and the Identification and Analysis of Clay Minerals*. Oxford University Press. Link to the second edition (1997). Bibliographic record for this edition.
- Drits, V. A.; Tchoubar, C., 1990. *X-Ray Diffraction by Disordered Lamellar Structures: Theory and Applications to Microdivided Silicates and Carbons*. Springer.

- Lanson, B.; Bouchet, A., 1995. *Identification des minéraux argileux par diffraction des rayons X : apport du traitement numérique*. Bulletin des Centres de Recherches Exploration-Production Elf Aquitaine, 19(1), 91–118. Full-text PDF hosted by ISTerre.
- Meunier, A., 2005. *Clays*. Springer.
- Poppe, L. J.; Paskevich, V. F.; Hathaway, J. C.; Blackwood, D. S., 2001. *Clay Mineral Identification Flow Diagram*. *A Laboratory Manual for X-Ray Powder Diffraction*, USGS Open-File Report 2001–41. Flow diagram in PDF format.

Full text in Portuguese · References and original document

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BAILEY, S. W. Structures of Layer Silicates. In: BRINDLEY, G. W.; BROWN, G. (ed.). Crystal Structures of Clay Minerals and their X-Ray Identification. London: Mineralogical Society, 1980. Mineralogical Society Monograph, v. 5. DOI 10.1180/mono-5.1.

BAILEY, S. W. (ed.). Hydrous Phyllosilicates (Exclusive of Micas). Washington, D.C.: Mineralogical Society of America, 1988. Reviews in Mineralogy, v. 19. Catálogo e sumário da editora.

BRINDLEY, G. W.; BROWN, G. (ed.). Crystal Structures of Clay Minerals and their X-Ray Identification. London: Mineralogical Society, 1980. Mineralogical Society Monograph, v. 5. DOI 10.1180/mono-5.

DRITS, V. A.; TCHOUBAR, C. X-Ray Diffraction by Disordered Lamellar Structures: Theory and Applications to Microdivided Silicates and Carbons. Berlin; Heidelberg: Springer-Verlag, 1990. DOI 10.1007/978-3-642-74802-8.

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Additional bibliography

References retained from the previous guide. Check the original publication before citing.

1. Basics of X-ray diffraction Stanjek, H., Häusler, W. Hyperfine Interactions, 2004
2. The routine mineralogical analysis of clay-bearing samples. Starkey, H.C., Blackmon, P.D., Hauff, P.L. US Geological Survey Bulletin, 1984
3. Quantitative analysis of clay and other minerals in sandstones by X-ray powder diffraction (XRPD) Hillier, S. Clay Mineral Cements in Sandstones, 2009
4. Quantitative analysis of clay and other minerals in sandstones by X-ray powder diffraction (XRPD) Hillier, S. International Association of Sedimentologists Special Publications, 2002
5. A comparative quantitative study of an Omani soil using X-ray diffraction technique Al-Rawas, A.A., Hago, A., Sutherland, H.H., (...), Al-Shihi, B. Geotechnical and Geological Engineering, 2001
6. Quantitative X-ray diffraction analysis of clay-bearing rocks from random preparations Środoń, J., Drits, V.A., McCarty, D.K., (...), Eberl, D.D. Clays and Clay Minerals, 2001
7. Quantitative X-ray powder diffraction analyses of clays using an orienting internal standard and pressed disks of bulk shale samples Cody, R.D., Thompson, G.L. Clays and Clay Minerals, 1976
8. Identification and characterization of clays using specific X-ray diffraction and computer modeling Hluchy, M. Ceramic Engineering and Science Proceedings, 2001
9. Quantitative analysis method of bulk mineralogy: Take the 9th Reynolds Cup samples for example Li, Y., Liu, Z. Marine Geology and Quaternary Geology, 2025
10. AUTOMATED FULL-PATTERN SUMMATION OF X-RAY POWDER DIFFRACTION DATA FOR HIGH-THROUGHPUT QUANTIFICATION OF CLAY-BEARING MIXTURES M. Butler, B., Hillier, S. Clays and Clay Minerals, 2021
11. New approaches to XRD profile modeling Delonas, C. Spectroscopy (Santa Monica), 2020
12. Machine-learning-based semi-quantitative analysis model of clay minerals in Onsan–Busan coastal sediments, South Korea Park, J.Y., Jo, Y., Koo, H.J., (...), Lee, K. Geosciences Journal, 2026
13. Quantitative mineralogical analysis of small samples of china clay using x ray diffractometry Salt, P.D. British Journal of Industrial Medicine, 1985
14. A Rapid Preparation Technique for X-Ray Diffraction Analysis of Clay Minerals in Weathered Rock Materials Mcalister, J.J., Smith, B.J. Microchemical Journal, 1995
15. X-ray diffraction techniques Whittig, L.D., Allardice, W.R. Methods of Soil Analysis, Part 1: Physical and Mineralogical Methods, 2018
16. Recognition of interstratified clays Mills, J.G., Zwarich, M.A. Clays and Clay Minerals, 1972
17. High-resolution X–ray diffraction datasets: Carbonates Amao, A.O., Al-Otaibi, B., Al-Ramadan, K. Data in Brief, 2022

18. The application of portable X-ray diffraction to quantitative mineralogical analysis of hydrothermal systems Burkett, D.A., Graham, I.T., Ward, C.R. Canadian Mineralogist, 2015
19. Structural analysis of interstratified illite-smectite by the rietveld method Wang, X., Wang, H. Crystals, 2021
20. Consistent modeling of the XRD patterns of mixed-layer phyllosilicates Plançon, A. Clays and Clay Minerals, 2004
21. Expert system for structural characterization of phyllosilicates: I. Description of the expert system Plancon, A., Drits, V.A. Clay Minerals, 1994
22. Mineralogical, geochemical characterization and application of an alluvial clay material Katte, V.Y., Ngapgue, F., Tsang, J.M. SN Applied Sciences, 2019
23. Distribution of clay minerals along a weathered pyroclastic profile, Hong Kong Duzgoren-Aydin, N.S., Aydin, A., Malpas, J. Catena, 2002
24. FREEZE DRYING: AN EFFECTIVE METHOD FOR PREPARING RANDOM X-RAY-DIFFRACTION MOUNTS FROM SMALL QUANTITIES OF CLAYS Maithel, S.A., Nick, K.E. Journal of Sedimentary Research, 2026
25. A Comparison Study of Synkinematic Illite Isolation, Quantitative X-ray Powder Diffraction, and K-Ar Dating for Direct Fault Gouge Analyses Zheng, Y., Haibing, L.I., Junjie, L.I., (...), Jialiang, S.I. Acta Geologica Sinica (English Edition), 2023
26. Use of X-ray transmission diffractometry for the study of clay-particle orientation of different water contents Inigo, A.C., Tessier, D., Pernes, M. Clays and Clay Minerals, 2000
27. Variation of preferred orientation in oriented clay mounts as a result of sample preparation and composition Dohrmann, R., Rüping, K.B., Kleber, M., (...), Jahn, R. Clays and Clay Minerals, 2009
28. Results of an inter-laboratory comparison of methods for quantitative clay analysis Ottner, F., Gier, S., Kuderna, M., Schwaighofer, B. Applied Clay Science, 2000
29. Utilizing smear mounts for X-ray diffraction as a fully quantitative approach in rapidly characterizing the mineralogy of shale gas reservoirs Munson, E.O., Chalmers, G.R.L., Bustin, R.M., Li, K. Journal of Unconventional Oil and Gas Resources, 2016
30. Microwave drying of clays for X-ray diffraction analysis. McCallister, D.L. Soil Science Society of America Journal, 1986
31. Fired pressed pellet as a sample preparation technique of choice for an energy dispersive X-ray fluorescence analysis of raw clays Mijatović, N., Vasić, M., Miličić, L., (...), Radojević, Z. Talanta, 2023
32. Sample Preparation for Energy-Dispersive X-Ray Fluorescence Analysis of Archaeological and Geochemical Silicic Samples: Selection of Pressed Powder Pellet or Loose Powder Specimens Ichikawa, S., Kitajima, N., Okubo, H., Kurisaki, T. X-Ray Spectrometry, 2026
33. Loose-powder technique for X-ray fluorescence analysis of ancient pottery using a small (100mg) powdered sample Ichikawa, S., Nakayama, K., Nakamura, T. X-Ray Spectrometry, 2012
34. X-ray diffraction of the five component clay shales Tangchawal, Sanga null, 1991

35. INTERSTRAT - an expert system to help identify interstratified clay minerals from powder XRD data: II. Testing the program Garvie, L.A.J. Clay Minerals, 1994

36. Interparticle diffraction: a new concept for interstratified clays. Nadeau, P.H., Wilson, M.J., McHardy, W.J., Tait, J.M. Clay Minerals, 1984

37. Interstratified clays as fundamental particles Nadeau, P.H., Wilson, M.J., McHardy, W.J., Tait, J.M. Science, 1984

38. Fourier transform methods for studying scattering from lamellar systems - I. A direct method for analysing interstratified mixtures MacEwan, D.M.C. Kolloid-Zeitschrift, 1956

39. Spectral identification (1.2-2.5 µm) and characterization of Paris Basin kaolinite/smectite clays using a field spectrometer Kruse, F.A., Thiry, M., Hauff, P.L. Physical measurements and signatures in remote sensing. Proc. 5th international colloquium, Courchevel, 1991. Vol. 1, 1991

40. Microanalysis of clay-based pigments in paintings by XRD techniques Hradil, D., Bezdička, P., Hradilová, J., Vašutová, V. Microchemical Journal, 2016

41. Energy dispersive X-ray diffraction and fluorescence portable system for cultural heritage applications Mendoza Cuevas, A., Bernardini, F., Gianoncelli, A., Tuniz, C. X-Ray Spectrometry, 2015